



PESIEREM1\_I NOX

SETS OF STAINLESS STEEL WEIGHTS IN M1 PRECISION CLASS




Set of STAINLESS STEEL weights suitable for use in legal metrology and industrial applications. The M1 precision class allows for the use of these weights for the check and calibration of scales and instruments in class III up to 10000e.

Our solutions for quality assurance and certainty of your measurements.

## MAIN FEATURES

- In compliance with the OIML - R111 recommendation;
- M1 precision class;
- in highly polished austenitic stainless steel;
- 7900 kg/m<sup>3</sup> density (except for fraction weights);
- initial verification stamp available upon request;
- Accredia (ea) calibration certificate upon request;
- wooden case;
- fitted with glove and/or tweezer (tweezer and dust brush with the WSAM1M05)



## LASER MARKING SERVICE UPON REQUEST:

- The laser marking service is useful to identify each sample weight or mass individually. More specifically:
    - The marking allows easier identification and traceability of each weight or mass;
    - It is perfectly smooth, does not remove material from the weight and does not affect the surface, thus avoiding dust deposit;
    - It is customisable, according to the requirements of the companies.
- In addition:
- It remains unchanged over time;
  - It can be performed on previously purchased masses and weights, of all accuracy classes, from 1 g to 20 kg;
  - It is perfect for pharmaceutical industries and laboratories.

## DETAIL 1



Example of optional laser marking

## VERSIONS

### WEIGHT SETS M1 CLASS

Codice	Numero pezzi	Pesi contenuti	Contenuto totale
WSAM1M05	12	mg: 1, 2, 2, 5, 10, 20, 20, 50, 100, 200, 200, 500	1,11g
WSAM1M50	20	mg: 1, 2, 2, 5, 10, 20, 20, 50, 100, 200, 200, 500 g: 1, 2, 2, 5, 10, 20, 20, 50	111,11g
WSAM1M100	21	mg: 1, 2, 2, 5, 10, 20, 20, 50, 100, 200, 200, 500 g: 1, 2, 2, 5, 10, 20, 20, 50, 100	211,11g
WSAM1M200	23	mg: 1, 2, 2, 5, 10, 20, 20, 50, 100, 200, 200, 500 g: 1, 2, 2, 5, 10, 20, 20, 50, 100, 200, 200	611,11g
WSAM1M500	24	mg: 1, 2, 2, 5, 10, 20, 20, 50, 100, 200, 200, 500 g: 1, 2, 2, 5, 10, 20, 20, 50, 100, 200, 200, 500	1111,11g
WSAM1MK1	25	mg: 1, 2, 2, 5, 10, 20, 20, 50, 100, 200, 200, 500 g: 1, 2, 2, 5, 10, 20, 20, 50, 100, 200, 200, 500, 1000	2111,11g
WSAM1MK2	27	mg: 1, 2, 2, 5, 10, 20, 20, 50, 100, 200, 200, 500 g: 1, 2, 2, 5, 10, 20, 20, 50, 100, 200, 200, 500, 1000, 2000, 2000	6111,11g
WSAM1MK5	28	mg: 1, 2, 2, 5, 10, 20, 20, 50, 100, 200, 200, 500 g: 1, 2, 2, 5, 10, 20, 20, 50, 100, 200, 200, 500, 1000, 2000, 2000, 5000	11111,11g
WSAM1MK10	29	mg: 1, 2, 2, 5, 10, 20, 20, 50, 100, 200, 200, 500 g: 1, 2, 2, 5, 10, 20, 20, 50, 100, 200, 200, 500, 1000, 2000, 2000, 5000, 10000	21111,11g
WSAM150	8	g: 1, 2, 2, 5, 10, 20, 20, 50	110g
WSAM1100	9	g: 1, 2, 2, 5, 10, 20, 20, 50, 100	210g
WSAM1200	11	g: 1, 2, 2, 5, 10, 20, 20, 50, 100, 200, 200	610g
WSAM1500	12	g: 1, 2, 2, 5, 10, 20, 20, 50, 100, 200, 200, 500	1110g
WSAM1K1	13	g: 1, 2, 2, 5, 10, 20, 20, 50, 100, 200, 200, 500, 1000	2110g
WSAM1K2	15	g: 1, 2, 2, 5, 10, 20, 20, 50, 100, 200, 200, 500, 1000, 2000, 2000	6110g
WSAM1K5	16	g: 1, 2, 2, 5, 10, 20, 20, 50, 100, 200, 200, 500, 1000, 2000, 2000, 5000	11110g
WSAM1K10	17	g: 1, 2, 2, 5, 10, 20, 20, 50, 100, 200, 200, 500, 1000, 2000, 2000, 5000, 10000	21110g



QUALITY ASSURANCE AND CERTAINTY OF YOUR MEASUREMENTS



**OFFICE**

Via Picasso, 18/20  
20025 Legnano (MI) - Italy



SALES AND TECHNICAL ASSISTANCE SERVICE

*All the specifications and images can be subject to variations and upgrades without prior notice.*